

Drying Oven



GZX(300°C) Series



BGZ(300°C) Series

Shanghai Boxun Medical Biological Instrument Corp
上海博迅医疗生物仪器股份有限公司

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com




www.sh-bx.com

Drying Oven >>

GZX (300°C) Series

Product Feature

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted at will;
- PMMA I control system, LCD display, equipped with temperature control, timing, and over temperature alarm functions;
- Preheating chamber design with power supply recovery function;
- Door with rotating two-stage locking structure.



Drying Oven >>

BGZ (300°C) Series

Product Feature

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted at will;
- PMMA II operating system, LCD display of various parameters, with temperature control, timing, programming, wind speed control, over temperature alarm and other functions;
- Door with rotating two-stage locking structure;
- Preheating chamber design, with power supply recovery function, appointment function;
- 6-segment adjustable wind speed, programmable, 10 segments with 100 cycles.



Stainless Steel Inner CaVity



Seal Ring



Doorknob



Control Panel

Model	GZX-9076MBE	GZX-9146MBE	GZX-9246MBE	GZX-9306MBE	GZX-9456MBE
Volume [L]	64	129	210	280	452
Temp. Range [°C]	RT+5~300	RT+5~300	RT+5~300	RT+5~300	RT+5~300
Temp. Resolution [°C]	0.1	0.1	0.1	0.1	0.1
Temp. Fluctuation at 150°C [°C]	±1	±1	±1	±1	±1
Temp. Uniformity at 150°C [%]	±2.5	±2.5	±2.5	±2.5	±2.5
Heating Rate to 150°C [°C/min]	>5	>5	>5	>5	>5
Power [W]	1550	2150	2750	2750	3500
Power Supply [V, Hz]	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%
Timer Range [min or h]	0~9999	0~9999	0~9999	0~9999	0~9999
Internal Dimensions W×D×H [mm]	410×350×450	570×400×570	630×530×630	720×530×745	600×580×1300
Housing Dimensions W×D×H [mm]	555×566×722	715×616×842	775×746×902	865×746×1017	774×845×1704
Number of Shelves (std./max.) [pc]	2 / 5	2 / 7	3 / 8	3 / 10	4 / 18
Load per Shelf [kg]	15	15	15	15	15
Net Weight [kg]	42	57	77	97	145
Gross Weight [kg]	60	75	96	118	174
Package Size [mm]	675*726*872	835*776*992	895*906*1052	985*906*1167	894*1005*1854

Model	BGZ-76	BGZ-146	BGZ-246	BGZ-306	BGZ-456	BGZ-656	BGZ-1006
Volume [L]	64	129	210	280	452	650	1000
Temp. Range [°C]	RT+5~300	RT+5~300	RT+5~300	RT+5~300	RT+5~300	RT+5~300	RT+5~300
Temp. Resolution [°C]	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Temp. Fluctuation at 150°C [°C]	±0.8	±0.8	±0.8	±0.8	±0.8	±0.8	±0.8
Temp. Uniformity at 150°C [%]	±2	±2	±2	±2	±2	±2	±2
Heating Rate to 150°C [°C/min]	>5	>5	>5	>5	>5	>5	>5
Power [W]	1550	2150	2750	2750	3500	4500	6600
Power Supply [V, Hz]	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	220±10%, 50±2%	380±10%, 50±2%
Timer Range [min or h]	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999
Reservation Range [min]	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999
Internal Dimensions W×D×H [mm]	410×350×450	570×400×570	630×530×630	720×530×745	600×580×1300	720×650×1400	1100×650×1400
Housing Dimensions W×D×H [mm]	555×566×722	715×616×842	775×746×902	865×746×1017	774×845×1704	865×866×1730	1220×850×1730
Number of Shelves (std./max.) [pc]	2 / 5	2 / 7	3 / 8	3 / 10	4 / 18	4 / 18	4 / 18
Load per Shelf [kg]	15	15	15	15	15	15	15
Net Weight [kg]	42	57	77	97	145	-	-
Gross Weight [kg]	60	75	96	118	174	-	-
Package Size [mm]	675*726*872	835*776*992	895*906*1052	985*906*1167	894*1005*1854	-	-